


# FORCEPS APPLICATION



# FORCEPS APPLICATION REFERENCE CHART

## UPPER FORCEPS

	INCISORS	CANINES	PREMOLARS	1st MOLARS	2nd MOLARS	3rd MOLARS	ROOTS
1	●	●					
10S				●	●	●	
18L				●	●		
18R				●	●		
53L				●	●		
53R				●	●		
62	●	●	●				
65	●						●
69							●
88L				●	●		
88R				●	●		
99C	●	●	●				
101			●				
150	●	●	●				●
150A	●	●	●				
150S*	●	●	●	●	●	●	●
210H					●		
210S					●		
286	●		●				●

(blotted broken teeth)

S = pediat etc

(narrow)

## LOWER FORCEPS

16				●	●		
17				●	●		
23				●	●		
62	●	●	●				
69							●
101			●				
151	●	●	●				●
151A	●	●	●				
151S*	●	●	●	●	●	●	●
203	●	●	●				●
222					●		

\*Pedodontic.

Note: Forceps Nos. 62 & 101 may be used for deciduous or permanent teeth.

For more detailed information on individual forceps, see your Hu-Friedy catalog.

Rt tip - rangeurs